







EDAN Edan Instruments, Inc.

- 3/F B, Nanshan Medical Equipments Park, Nanhai Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China
- Tel: +86-755-26898326 Fax: +86-755-26898330 www.edan.com.cn Email: info@edan.com.cn
- All rights reserved. Features and specifications are subject to change without notice







SD3 VASCULAR



SD3 LITE







Configurations

Model	SD3 VASCULAR	SD3 LITE	SD3	SD3 PLUS	SD3 PRO
FHR Display	×	×	$\sqrt{}$	$\sqrt{}$	$\sqrt{}$
Charger Stand	×	×	×	0	0
Record & Play	×	×	×	×	$\sqrt{}$
AA Batteries	0	0	0	×	×
Li-ion Battery	×	×	×	V	V